

# 5Ws\* OMNICHROMA *Flow* BULK

This composite resin uses smart chromatic technology to shade match any tooth color from A1 to D4. [compiled by Kristin Hohman]

Information provided by Tokuyama Dental America.

## THE 5Ws\*

### WHAT

OMNICHROMA *Flow* BULK composite resin

### WHO

Tokuyama Dental America

877-378-3548 | [omnichromaflowbulk.com](http://omnichromaflowbulk.com)

CIRCLE RS #22

### WHERE

All dental practices that perform restorative care looking to save time, reduce cost, and minimize inventory.

### WHEN

Indicated for direct anterior and posterior restorations, OMNICHROMA *Flow* BULK can also be used as a cavity base or liner and to block out cavity undercuts before fabricating indirect restorations. It can also be used in the repair of porcelain or composite.

### WHY

OMNICHROMA *Flow* BULK is a low-viscosity, light-cured, radiopaque composite resin for use in direct restorations. It is a single-shade matching material. OMNICHROMA *Flow* BULK is a flowable material that can be placed in 3.5-mm increments and is available in syringe form.



## THE HOW

Employing smart chromatic technology and uniformly sized spherical fillers, OMNICHROMA *Flow* BULK composite generates a red to yellow color range to esthetically shade match any tooth color from A1 to D4. The color generation ensures there is no loss of chroma or drop in value. The composite's fillers are formulated to generate color, reflecting the shade of the surrounding tooth. OMNICHROMA *Flow* BULK also features low polymerization shrinkage, high compressive strength, excellent cavity adaptation, and high polishability.